



CAUTION Investigational Device.  
Limited by Federal (or United States)  
law to investigational use.

1609 Shoal Creek Blvd., Ste. 303  
Austin, TX 78701 USA

Patient Name:

Patient ID:

### For Medical Questions, Contact Your Physician

Surgical Center:

Phone:

This person is implanted with an: **EndoStim LES Stimulation System – LES Stimulator and Stimulating Lead**

Stimulator SN:

Lead SN:

PN 6600 Rev. C 02/19/2018

The following treatments are contraindicated (not recommended) when implanted with the EndoStim LES Stimulation System:

- Medical Diathermy
- TENS Exposure in the Abdominal Region

**MR Conditional.** Contact your EndoStim Physician before scheduling an Magnetic Resonance Imaging (MRI) procedure.



In an emergency, placement of a magnet over the device for at least two (2) seconds terminates therapy delivery for 24-48 hours.